



X-Ray Protective Eyewear – Portfolio IPL 2016/2

Valid December 1st, 2015

MAVIG

General Notice

All sales transactions are done in EURO. Unless otherwise stated, prices of goods includes packing, and are FCA MAVIG, Munich (Incoterms 2010). In order to protect the interest of the purchaser, we recommend to obtain insurance to cover against loss, damage, and or theft of goods in transit.

Returned merchandise will not be accepted without prior written notice. Any merchandise received by MAVIG without written authorization shall be refused and returned to sender.

Special orders are sold on a non-refundable basis. Aprons ordered without colour specification will be supplied in the standard colour listed for each model. All people authorized to use MAVIG products should be trained in their correct use and storage, as well as the dangers involved in exposure to radiation. The proper selection and use of MAVIG products are the sole responsibility of the purchaser. The purchase must only be made by a person competent in determining the requirements of X-Ray protection and the particular applications in which they will be used.

No warranty or guarantee, either expressed or implied, is extended to cover personal injury or damage which may result due to improper use, storage, handling, or unauthorized use of MAVIG products. This warranty is in lieu of all warranties expressed or implied, including warranties of merchantability and fitness for purpose.




We refer to the general terms of business of MAVIG GmbH.

Right to technical changes for the products or content of the catalogue is reserved.




X-Ray Protective Eyewear 2016/2 - Product Overview

Page	Model	Type	Design / Colour	incl. Accessories			Prescription & Coating			
				Case	Cleaning Cloth	Band	Single	Bifocal	Progressive	Anti-Reflection
4	BR126	Glasses Glasses	Cocoa with Light Blue Silver Bronze with Orange	• •	• •	• •	• •	- -	• •	• •
5	BR115	Glasses	Black	•	-	•	-	-	-	-
5	BR118	Glasses	Blue, Taupe, Black	•	-	•	•	•	-	-
5	BR119	Glasses	Silver, Black	•	-	•	•	•	-	-
5	BR178	Glasses	White-transparent	•	-	•	-	-	-	-
6	BR310	Glasses	Tortoise	•	-	•	•	•	-	-
6	BR321	Glasses	Bronze	•	-	•	•	•	-	-
6	BR322	Glasses	Gold	•	-	•	•	•	-	-
6	BR331	Glasses	Blue-transparent	•	-	•	•	•	-	-
7	BRV500	Visor	Continuous protection from forehead to chin	-	-	-	-	-	-	-
7	BRV501	Visor	Protection down to the cheek bones with a cut-out for nose	-	-	-	-	-	-	-

X-Ray Protective Eyewear 2016/2

BR126	Article-No.	Size	Pb Glas	Prescription	
 <p>Design: Cocoa/Light Blue</p> <p>Weight: ca. 70 g (Pb 0.50 mm*) ca. 80 g (Pb 0.75 mm*)</p>  	with Anti-Reflection				
	BR126SM5C-arc BR126SM7C-arc	small-medium small-medium	0.50 mm 0.75 mm	plano plano	
	BR126SM5C-arc-presc-S BR126SM5C-arc-presc-P	small-medium small-medium	0.50 mm 0.50 mm	Single Progressive	
	BR126SM7C-arc-presc-S BR126SM7C-arc-presc-P	small-medium small-medium	0.75 mm 0.75 mm	Single Progressive	
	BR126L5C-arc BR126L7C-arc	large large	0.50 mm 0.75 mm	plano plano	
	BR126L5C-arc-presc-S BR126L5C-arc-presc-P	large large	0.50 mm 0.50 mm	Single Progressive	
	BR126L7C-arc-presc-S BR126L7C-arc-presc-P	large large	0.75 mm 0.75 mm	Single Progressive	

* The model BR126 is fitted with special side shields in Pb 0.50 mm.

BR126	Article-No.	Size	Pb Glas	Prescription	
 <p>Design: Silver Bronze/Orange</p> <p>Weight: ca. 70 g (Pb 0.50 mm*) ca. 80 g (Pb 0.75 mm*)</p>  	with Anti-Reflection				
	BR126SM5S-arc BR126SM7S-arc	small-medium small-medium	0.50 mm 0.75 mm	plano plano	
	BR126SM5S-arc-presc-S BR126SM5S-arc-presc-P	small-medium small-medium	0.50 mm 0.50 mm	Single Progressive	
	BR126SM7S-arc-presc-S BR126SM7S-arc-presc-P	small-medium small-medium	0.75 mm 0.75 mm	Single Progressive	
	BR126L5S-arc BR126L7S-arc	large large	0.50 mm 0.75 mm	plano plano	
	BR126L5S-arc-presc-S BR126L5S-arc-presc-P	large large	0.50 mm 0.50 mm	Single Progressive	
	BR126L7S-arc-presc-S BR126L7S-arc-presc-P	large large	0.75 mm 0.75 mm	Single Progressive	


* The model BR126 is fitted with special side shields in Pb 0.50 mm.


ATTENTION


For orders including corrective lenses, please indicate the following values from a current prescription:



- Diopter values (Sph)
- If applicable, cylindrical values (Cyl) and axis values (Axis)
- Pupillary distance, separately for your left and right eye (PD)
- Addition values (ADD) for bifocal or progressive lenses

The prescription should not be older than 6 months. Only after receiving all prescription values needed, can we confirm if the desired corrective lenses are possible.


Article-No.	Size	Pb Glas	Prescription	(without antireflection) BR115
BR115-Black	universal	0.75 mm	Plano	 Design: Black Weight: ca. 110 g (Pb 0.75 mm*)
Corrective lenses not available.				

Article-No.	Size	Pb Glas	Prescription	(without antireflection) BR118
BR118-Blue	universal	0.75 mm	Plano	 <p>Designs: Blue, Taupe, Black Weight: ca. 60 g (Pb 0.75 mm)</p>
BR118-Taupe	universal	0.75 mm	Plano	
BR118-Black	universal	0.75 mm	Plano	
BR118-99*	universal	0.75 mm	Single	
BR118-99*	universal	0.75 mm	Bifocal	
* Please indicate the desired colour design and current correction values when placing your order - more information under "ATTENTION".				

Article-No.	Size	Pb Glas	Prescription	(without antireflection) BR119
BR119-Silver	universal	0.75 mm	Plano	 <p>Designs: Silver, Black Weight: ca. 65 g (Pb 0.75 mm)</p>
BR119-Black	universal	0.75 mm	Plano	
BR119-99*	universal	0.75 mm	Single	
BR119-99*	universal	0.75 mm	Bifocal	
* Please indicate the desired colour design and current correction values when placing your order - more information under "ATTENTION".				

Article-No.	Size	Pb Glas	Prescription	(without antireflection) BR178
BR178*	universal	0.75 mm	Plano	  Design: White-transparent Weight: ca. 100 g (Pb 0.75 mm)
* Please note that the BR178 is not CE/PPE certified.				
Corrective lenses not available.				


X-Ray Protective Eyewear 2016/2

BR310 (without antireflection)	Article-No.	Size	Pb Glas	Prescription	
	BR310	universal	0.75 mm	Plano	
	BR310-99*	universal	0.75 mm	Single	
	BR310-99*	universal	0.75 mm	Bifocal	
<p>* Please indicate current correction values when placing your order - more information under "ATTENTION".</p>					

Design: Tortoise

Weight: ca. 90 g (Pb 0.75 mm*)


* Side protection Pb 0.75 mm

BR321 (without antireflection)	Article-No.	Size	Pb Glas	Prescription	
	BR321	universal	0.75 mm	Plano	
	BR321-99*	universal	0.75 mm	Single	
	BR321-99*	universal	0.75 mm	Bifocal	
<p>* Please indicate current correction values when placing your order - more information under "ATTENTION".</p>					

Design: Bronze

Weight: ca. 65 g (Pb 0.75 mm*)


* Side protection Pb 0.50 mm

BR322 (without antireflection)	Article-No.	Size	Pb Glas	Prescription	
	BR322	universal	0.75 mm	Plano	
	BR322-99*	universal	0.75 mm	Single	
	BR322-99*	universal	0.75 mm	Bifocal	
<p>* Please indicate current correction values when placing your order - more information under "ATTENTION".</p>					

Design: Gold

Weight: ca. 70 g (Pb 0.75 mm*)

* Side protection Pb 0.50 mm

BR331 (without antireflection)	Article-No.	Size	Pb Glas	Prescription	
	BR331	universal	0.75 mm	Plano	
	BR331-99*	universal	0.75 mm	Single	
	BR331-99*	universal	0.75 mm	Bifocal	
<p>* Please indicate current correction values when placing your order - more information under "ATTENTION".</p>					

Design: Blue-transparent

Weight: ca. 85 g (Pb 0.75 mm*)

* Side protection Pb 0.50 mm

Article-No.	Size	Pb Glas	Prescription	
BRV500	universal	0.10 mm	Plano	

Corrective lenses not available.

Shield made of lead acrylic,
without antireflection.



Design: Continuous protection from
forehead to chin

Weight: ca. 425/525 g (without/with pads)

Article-No.	Size	Pb Glas	Prescription	
BRV501	universal	0.10 mm	Plano	

Corrective lenses not available.

Shield made of lead acrylic,
without antireflection.



Design: Protection down to the cheek bones
with a cut-out for the nose

Weight: ca. 405/505 g (without/with pads)

ATTENTION

*For orders including corrective lenses, please indicate
the following values from a current prescription:*

- Diopter values (Sph)
- If applicable, cylindrical values (Cyl) and axis values (Axis)
- Pupillary distance, separately for your left and right eye (PD)
- Addition values (ADD) for bifocal or progressive lenses

The prescription should not be older than 6 months. Only after receiving all prescription values needed,
can we confirm if the desired corrective lenses are possible.



MAVIG GmbH Headquarters

PO Box 82 03 62
81803 Munich
Germany

Stahlgruberring 5
81829 Munich
Germany

Phone: +49 (0) 89 420 96 0
Fax: +49 (0) 89 420 96 200
e-Mail: info@mavig.com

WORLDWIDE MAVIG OFFICES

Nordic & Baltic Countries MAVIG Nordic

Stockholm
Sweden

Phone: +46 (0) 722 25 25 68
e-Mail: larsson@mavig.com

Benelux, GB, Ireland MAVIG B.V.

Dunne Bierkade 28
2512 BD Den Haag
Netherlands

Phone: +31 (0) 70 345 99 85
Mobile: +31 (0) 61 595 43 48
e-Mail: simmonds@mavig.nl

France MAVIG France SAS

66, Ave. des Champs Elysées
F-75008 Paris
France

Phone: +33 (0)1 30 59 46 23
Fax: +33 (0)1 30 59 46 23
e-Mail: info@mavig.fr

China MAVIG Healthcare Beijing Co., Ltd.

Tower 1, Block C - Room No. 2607
Wangjing SOHO
Futong East Street No. 1
Wangjing Area, Chaoyang District
Beijing 100102, China

Phone: +86 - 10 - 46 79 62 88
Fax: +86 - 10 - 46 79 62 88
Mobile: +86 - 13 70 122 85 86
e-Mail: zhang@mavig.com

REPRESENTATION

USA & Canada Ti-Ba Enterprises, Inc.

25 Hytec Circle
Rochester, NY 14606
USA

Phone: +1 (1) 585 247 1212
Fax: +1 (1) 585 247 1395
(international)
1 800 836 8422
(in the USA)
e-Mail: mavig@ti-ba.com

Middle East & Northern Africa Mena Medical Development

Ashrafieh - Adlieh Square - Alfaras St.
- Alboustany Building, 5th floor
Beirut
Liban

Phone: +961 14 23 499
Fax: +961 14 26 499
e-Mail: mavig@mena-md.com

www.mavig.com